

## Colony PCR

1. After transformations, colonies are growing on the plate next day.
2. Prepare one new plate with specific antibiotics, dry and warm up to 37 degree.
3. Prepare PCR master mix (250 ul) and aliquot 25 ul to each of 10 PCR tubes

10 X buffer	25 ul
dNTP (25 mM)	2.5 ul
Primer F (20 uM)	2.5 ul
Primer R (20 uM)	2.5 ul
Taq	1.0 ul
H2O	216.5 ul
4. Pipette 0.5 ul positive control (DNA) to PCR tube 1
5. Pick up 9 colonies from the old plate by pipette tips, dip into PCR tube and strip on the new plate with proper mark (see the pictures).
6. Place the new plate in 37 degree incubator for cells to grow
7. Run standard PCR protocols
8. Run samples on agarose gel
9. Select several colonies which show desire bands and culture in the liquid medium at 37 degree overnight.

